

Aluminium Over Toilet Frame

In some homes the toilet may be positioned too low for those who have difficulty bending or lowering themselves, particularly those who may be requiring or recovering from a hip replacement. Raising the height of the toilet seat makes toileting much safer & easier. Over Toilet Frames with a seat & lid, have arms for added stability, support & grip, while the legs are height adjustable.

The AusCare Aluminium over toilet frame assists transferring on & off the toilet. Constructed from sturdy no rust aluminium, this frame includes a lifetime warranty against rust.

The plastic seat & lid clip off to facilitate hygienic & easy cleaning. The anti-slip rubber tips offer stability & steadiness even on wet surfaces.

This aluminium over toilet frame is one of the very few on the market to be 100% Australian Standards compliant approved.

Easy push button height adjustment legs. The uncoated internal sliding legs provide smooth & effortless adjustment whilst avoiding any chance of water leakage into the legs.

The reinforced frame uses 25mm tubing, (larger than other standard frames), enabling users to gain a firmer grip whilst making transfers.

Specifications

Weight with lid	3.4kg
Weight without lid	2.7kg
Seat Width (Between Arms)	450mm
Overall Width (Footprint)	550mm
Overall Depth (Footprint)	550mm
Seat Height Range	440 - 590mm
Armrest Height (Seat to Top of Handle)	180mm
Australian Standards Tested	100kg
Manufacturer's Working Load	150kg
Height Clearance around Toilet	480mm
Height Clearance over Toilet	400 - 550mm



Warranty

A+ Lifetime Product Warranty against rust and corrosion. 12 month warranty on faulty workmanship and defective parts

Code	Description
NOV-AC02A	Over Toilet Frame
NOV-BAT68662	Bowl & Lid
NOV-8560-B01	Splashguard
NOV-N8500CA17	Toilet Seat
NOV-N8500CA18	Toilet Lid
NOV-BAT70373	Rubber Stoppers
NOV-BAT72012	Rubber Seat & Lid Stoppers