



2023
Product
Catalogue



Selecting the right mobility aids to meet your needs is an important decision. Many people use assistive devices to increase or maintain independence. Wheelchairs, rollators, walkers, crutches and walking sticks empower you to go further and fully experience the world around you — on your terms.

For over half a century, Days® has been a recognised global leader in empowering personal mobility. Fully understanding and embracing an individual's mobility needs, Days enhances and advances lives at every turn. When it comes to safety, ease of use, affordable quality, customer support and breadth of product choices, the world turns to Days.

When you select a Days mobility aid, you can expect to experience the difference. Days represents:

Boundlessness

We share our passion for removing life's boundaries, offering clinicians and individuals a broad range of mobility solutions. With a global reach, Days is always nearby.

Moving

Days is trust that moves you. We strive to inspire excitement about the possibilities of everything life has to offer. We believe your physical needs should not be a limitation, and we are here to help you propel your life — wherever you wish to go.

Knowledgeable

Days products and services support you by being built on a solid foundation of quality materials and heritage. With a firm understanding of the mobility industry, we are made to exceed expectations. You can rest assured with Days.

Compassionate

With every design, Days considers the unique needs and demands of the user. Days knows the products it creates will move both the body and spirit of an individual.



Wheelchairs 4

Rollators 14

Walking Frames 27

Crutches 29

Walking sticks 32

Go further with Days

Founded in 1968 in Bridgend, Wales, Days has led the industry in mobility for over 50 years. With a broad range of wheelchairs, rollators, walkers, crutches and walking sticks, Days is the brand to which the world turns to advance their lives. Our products are easy to use and ready when you need them while offering affordable quality. With a steadfast focus on safety, Days offers products that customers trust for their mobility needs.

The Days foundation rests on four central pillars: dedication, empowerment, understanding and enabling freedom.

Dedication

With a steadfast focus on product safety and ease of use, customers and clinicians can trust — and rest assured — knowing our products will take care of their mobility needs. With over a half-century of expertise, Days has a time-honoured reputation that continues to propel customers and their lives forward each day.

Empowerment

Days goes the distance with a quality, broad range of mobility aids that fit each unique individual for life's ever-changing path. Days helps empower clinicians, distributors and individuals to find the ideal solution for daily freedom of movement and motion.

Understanding

Days has developed an array of mobility support with every diverse need in mind. With extensive knowledge of the mobility industry, Days applies its empathy and understanding to produce quality, affordable products that support motion.

Enabling freedom

Navigating through life's twists and turns, Days helps keep the momentum going with a determination to uplift, support and answer each user's needs.



Wheelchairs

You're in the driver's seat

The World Health Organization recognizes that wheelchairs are essential for mobility, dignity and general wellbeing¹. Whether self- or attendant-propelled, your high-quality chair should enhance your mobility and safety while being easy to use and affordable. With Days, you will also have access to customer support when you need it most.

All Days Wheelchairs include the following features:

- Puncture-proof tyres to help ensure trouble-free mobility, low maintenance and a smooth ride over most surfaces
- Moulded armrests for extra comfort
- Folding for easy storage and transportation
- Wipe-clean padded seat fabric



References: 1. Strengthening access to quality wheelchair services. World Health Organization. <https://www.who.int/teams/health-product-policy-and-standards/assistive-and-medical-technology/assistive-technology/wheelchair-services>. Accessed May 27, 2022.

Escape Aluminium Wheelchairs

Attendant or Self-propelled

Ease of transportation and storage

Escape Wheelchairs are lightweight wheelchairs suitable for indoors and outdoors. Experience peace of mind and comfortability with puncture-proof tyres, a full-length armrest, padded upholstery and detachable, swing-away footrests. You will have the flexibility to get around no matter where you go.

Key features

- Lightweight, foldable frame with foldable backrest
- Swing-away and detachable footplates
- Padded, flame-retardant upholstery
- Includes a cable braking system
- Seatbelt included
- AS/NZS 3695.1:2011
- ARTG 200209

Specifications

Standard

Seat width	46 cm (18")
Seat depth	41 cm (16")
Seat height (front)	51.5 cm (20¼")
Seat height (rear)	49.5 cm (19½")
Width between arms	41 cm (16")
Backrest height above seat	40 cm (15¾")
Footrest length	39 - 49 cm (15 - 19")
Armrest height above seat	23 cm (9")

Attendant



Self-Propelled

External width (open)	53 cm (21")	64.5 cm (25")
External width (closed)	23 cm (9")	31 cm (12¼")
External depth (w/o footrests)	75 cm (29½")	80 cm (31½")
External depth (w/ footrests)	99.5 cm (39")	106 cm (41½")
Chair weight	12 kg (26 lbs)	15 kg (34 lbs)
Chair weight (w/o footrests)	10.5 kg (23 lbs)	14 kg (31 lbs)

Colour	Self-propelled	Attendant
Silver Blue	PAT-091566249	PAT-091171727
Ruby Red	PAT-091566256	PAT-091556893

<p>115 kg 18 st 253½ lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>2 YEAR Main frame warranty</p>	<p>Self-propelled</p>	<p>Attendant-propelled</p>
--	------------------------------------	---------------------------------------	-----------------------	----------------------------

Available colours



Wheelchairs

Swift Wheelchair

Attendant propelled

Our best selling wheelchair

Swift Wheelchairs are versatile, easy-to-use chairs. Lightweight with a robust and high-quality aluminium frame, Swift can be quickly disassembled without tools, making transportation and storage simple. Swift is compact, manoeuvrable and its sleek durable design makes it a great choice for rental pools and hospital use.

Key features

- Designed for indoor or outdoor occasional use
- Half-folding back for compact storage and convenient reaching
- Attendant-propelled version includes a cable braking system
- AS/NZS 3695.1:2011
- ARTG 200209

Specifications

Seat width	46 cm (18")
Seat depth	42 cm (16.5")
Seat height (front)	51 cm (20")
Seat height (rear)	48.5 cm (19")
Backrest height above seat	40 cm (15¾")
Footrest length	44 - 54 cm (17¼ - 21¼")
Armrest height above seat	23 cm (9")
External width (open)	57 cm (22½"), 62 cm (24½")
External width (closed)	24 cm (9½")
External depth (w/o footrests)	75 cm (29½")
External depth (w/ footrests)	107 cm (42")
Chair weight	11.5 kg (25.3 lbs)
Chair weight (w/o footrests)	10.2 kg (22.4 lbs)



Attendant propelled	Seat width
ANI-DAYWHE338T18	46 cm (18")
Optional leg rests	
ANI-DAYWHE338EL1-L	46 cm Elevating Legrest Standard - Left
ANI-DAYWHE338EL1-R	46 cm Elevating Legrest Standard - Right

115 kg 18 st 253½ lbs Maximum user weight	1 YEAR Product warranty	2 YEAR Main frame warranty	Attendant-propelled
---	----------------------------	-------------------------------	---------------------

Available colour

Cool Silver

Swift Wheelchair

Self-propelled

Our best selling wheelchair

Swift Wheelchairs are versatile, easy-to-use chairs. Lightweight with a robust and high-quality aluminium frame, Swift can be quickly disassembled without tools, making transportation and storage simple. Swift is compact, manoeuvrable and its sleek durable design makes it a great choice for rental pools and hospital use.

Key features

- Designed for indoor or outdoor occasional use
- Half-folding back for compact storage and convenient reaching
- Available in 3 seat widths
- AS/NZS 3695.1:2011
- ARTG 200209

Specifications

Seat width	35.5 cm (14"), 41 cm (16"), 46 cm (18")
Seat depth	42 cm (16.5")
Seat height (front)	51 cm (20")
Seat height (rear)	48.5 cm (19")
Backrest height above seat	40 cm (15¾")
Footrest length	44 - 54 cm (17¼ - 21¼")
Armrest height above seat	23 cm (9")
External width (open)	62 cm (24"), 67 cm (26")
External width (closed)	28 cm (11")
External depth (w/o footrests)	83 cm (32¾")
External depth (w/ footrests)	113 cm (44½")
Chair weight	12.5 kg (27.5 lbs)
Chair weight (w/o footrests)	11.2 kg (24.7 lbs)



Cable Handbrake System



Quick-release wheels

<p>115 kg 18 st 253½ lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>2 YEAR Main frame warranty</p>	<p>Self-propelled</p>	<p>3 Seat width options</p>
--	------------------------------------	---------------------------------------	-----------------------	---------------------------------

Available colour

Cool
Silver

Self-propelled	Seat width
ANI-DAYWHE338SP14	35.5 cm (14")
ANI-DAYWHE338SP16	41 cm (16")
ANI-DAYWHE338SP18	46 cm (18")
Optional leg rests	
ANI-DAYWHE338EL1-L	46cm Elevating Legrest Standard - Left
ANI-DAYWHE338EL1-R	46cm Elevating Legrest Standard - Right
ANI-DAYWHE338EL2	Paediatric Elevating Legrest - Pair
ANI-AAST-920000LH	Anti Tipper Wheels - Left
ANI-AAST-920000RH	Anti Tipper Wheels - Right
ANI-SP-DAYSW110002	Anti Tipper Rubber Stop

Wheelchairs

Premium Link Wheelchair

Self-propelled

Ideal for part- or full-time use

Height adjustable, lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Key features

- Light weight aluminium - 17.2kg
- Quick release wheels
- Solid tyres
- 16", 18", 20" and 22" seat widths available
- Swing in and out legrests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Drum brakes on the 20" and 22"
- Anti-tipper comes as standard inclusion
- ISO 7176-19
- ARTG 174573



Specifications

Seat width	41 cm (16"), 46 cm (18"), 51 cm (20"), 56 cm (22")
Seat depth	43 cm (17")
External width (open)	62 cm (24½")
External height (open)	98 cm (38½")
Chair weight	17.5kg (38.6 lbs)
Chair weight (w/o wheels)	14 kg (30.9 lbs)

Self-propelled	Seat width
ANI-DAYWHECWRI16	41 cm (16")
ANI-DAYWHECWRI18	46 cm (18")
ANI-DAYWHECWRI20D	51 cm (20")
ANI-DAYWHECWRI22D	56 cm (22")
Optional leg rests	
ANI-SP-DAYLINK0054	Amputee Stump Supports - Left
ANI-SP-DAYLINK0055	Amputee Stump Supports - Right
ANI-SP-DAYLINK0050	41-43cm Elevating Legrest Assembly - Left
ANI-SP-DAYLINK0052	41-43cm Elevating Legrest Assembly - Right
ANI-SP-DAYLINK0051	46-56cm Elevating Legrest Assembly - Left
ANI-SP-DAYLINK0053	46-56cm Elevating Legrest Assembly - Left

<p>16" & 18"</p> <p>130 kg</p> <p>20½ st 286¾ lbs</p> <p>Maximum user weight</p>	<p>20" & 22"</p> <p>150 kg</p> <p>23½ st 330½ lbs</p> <p>Maximum user weight</p>	<p>1 YEAR</p> <p>Product warranty</p>	<p>2 YEAR</p> <p>Main frame warranty</p>	<p>Self-propelled</p>	<p>4</p> <p>Seat width options</p>	<p>Crash-tested</p>
--	--	---------------------------------------	--	-----------------------	------------------------------------	---------------------

Available colour



Swift Lite Wheelchair

Attendant propelled

Push down armrests

Ultra Lightweight Wheelchair ideal for indoor and occasional use. The Swift Lite is designed with the carer in mind, the compact lightweight frame makes it easy to transfer in and out of a car. The push button height adjustable armrests make it versatile allowing you to push under a table then quickly adjust for comfort. Attendant propelled, in one size.

Key features

- Ultra lightweight frame 8kg
- Folding backrest
- Rear wheel locks
- 17" seat width
- Swing away, removable leg rests
- Height adjustable, flip up footplate
- ARTG 200209

Specifications

Seat width (between arms)	43 cm (17")
Seat depth	40 cm (15¾")
Seat height from floor	46 cm (18¼")
Backrest height above seat	40 cm (15¾")
Armrest height above seat	25.4 cm (10")
Legrest height range	43 - 53 cm (15¾ - 20¾")
External width	55 cm (21¾")
External depth (w footrests)	85 cm (33½")
Solid wheels	20.3 cm (8") rear, 15.2cm (6") front
Chair weight	9.5 kg (21 lbs)




Armrests adjust to fit under table

Colour	
Blue	ANI-ALTS1703-15

<p>115 kg 18 st 253½ lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>2 YEAR Main frame warranty</p>	<p>Attendant-propelled</p>
--	------------------------------------	---------------------------------------	----------------------------

Available colours



Wheelchairs

Whirl Steel Wheelchair

Attendant or Self-propelled

Robust and durable

Whirl Wheelchairs are simple yet versatile and ideal for part or full-time use. Whirl has a detachable, half-folding backrest for compact storage, and its narrow frame and compact dimensions ensure excellent manoeuvrability indoors and in confined spaces.

Key features

- Removable desk armrests to enable closer access to tables and worktops
- Solid tyres
- Swing away removable leg rests
- Self-propelled version features removable quick-release rear wheels
- ISO 7176-19
- ARTG 200209

Specifications

Seat width	45 cm (17¾")
Seat depth	42 cm (16½")
Seat height (front)	51 cm (20")
Seat height (rear)	48 cm (19")
Backrest height above seat	43 cm (17")
Footrest length	43.5 - 50 cm (17¼ - 19½")
Armrest height above seat	23 cm (9")

Self-propelled


External width (open)	Seat width +20 cm
External width (closed)	29 cm (11½")
External depth (w/o footrests)	80 cm (31½")
External depth (w/ footrests)	105 cm (41¼")
Chair weight	17 kg (37.4 lbs)
Chair weight (w/o footrests)	15 kg (33 lbs)

Attendant


Seat width +16 cm
25 cm (9¾")
70 cm (27½")
95 cm (37½")
15.4 kg (33.9 lbs)
13.4 kg (29.5 lbs)



Removable armrests for easier side transfer



Height-adjustable swing-away footrests



Easy to dismantle for storage or transport



Optional elevating leg rests



Available colour



Self-propelled	Attendant	Seat width
KUN-WHIRL45SP	KUN-WHIRL45TR	45 cm (17¾")
Optional leg rests		
KUN-SP-DAYWHIRL920004	43 - 45cm Elevating Legrest - Right	
KUN-SP-DAYWHIRL920002	43 - 45cm Elevating Legrest - Left	
KUN-SP-DAYWHIRL910001	Anti Tipper Wheels Self Propelled - Each	
KUN-WHIRL910000	Anti Tipper Wheels Self Propelled - Pair	

Days Tilt 'n' Space Wheelchair

Self-propelled

Comfort and function

This chair will suit most situations and is a must to consider for long patients requiring a high level of comfort and positioning options.

Key features

- Gas assisted tilt-in-space and recline options for optimum postural positioning
- Removable adjustable height arms, which feature adjustable arm pad length
- Removable swing away elevating leg rests and calf pads
- Height adjustable padded contoured backrest
- Removable padded headrest with adjustable side supports
- Alloy construction with water resistant breathable upholstery
- Tray comes as standard inclusion
- AS/NZS 3695.1:2011
- ARTG 174573

Specifications

Seat width	44 cm (17¼") : 49 cm (19¼")
Seat height (front)	51 cm (20")
Armrest height above seat	21 - 24 cm (8¼ - 9½")
Backrest height above seat	53.5 - 63.5 cm (21 - 25")
Seat board to footplate	42.5 - 61 cm (16¾ - 24")
Rear wheel size	24" x W28 mm
Front castor size	8" x W50 mm
Seat adjust angle	0 ~ 16 deg
Backrest adjust angle	90 ~ 135 deg
Footrest adjust angle	70 deg ~ horizontal
Net Weight	47 kg - 44 cm : 50kg - 49 cm



Tray included as standard



<p>136 kg 21 st 300 lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>Self-propelled</p>	<p>2 Seat width options</p>
---	------------------------------------	-----------------------	---------------------------------

Available colour



	Self-propelled
44 cm (17¼")	ANI-RC440
49 cm (19¼")	ANI-RC490
Replacment Tray	ANI-JM325

Wheelchairs

Transit Shopper Wheelchairs

Attendant-propelled



Features	Shopper 8 Transit	Shopper 12 Transit	Shopper 12 Transit
Folding Backrest	✓	✓	✓
Rear Wheel Locks	✓	✓	✓
Padded Armrests	✓	✓	✓
Safety Lap Belt	✓	✓	✓
Attendant Handbrakes	x	✓	✓
Height Adjustable Flip-Up Footplates	✓	✓	✓
Swing-Away Leg Rests	✓ Removable	✓ Removable	✓ Removable
Warranty	One Year Limited	One Year Limited	One Year Limited
ARTG	174474	174474	174474
Specifications			
Weight of Chair - Frame Only	9kg	11kg	13.5kg
Weight of Chair - with Leg Rests	10.5kg	12.5kg	15kg
Overall Seat Width	432mm	482.5mm	640mm
Seat Wide (between arms)	381mm	432mm	560mm
Seat Depth	407mm	407mm	407mm
Backrest Height	432mm	457mm	457mm
Leg Rest Height Range	356mm - 483mm	356 - 483mm	356 - 483mm
Overall Width (footprint)	508mm	559mm	680mm
Overall Depth (footprint)	584mm	635mm	940mm
Armrest Height (seat to top of handle)	254mm	254mm	254mm
Wheels (solid castors) - Front	203mm x 25.4mm	203mm x 25.4mm	203mm x 25.4mm
Wheels - Rear	203mm x 25.4mm	305mm x 57mm	505mm x 57mm
Colour	Black	Red	Blue
Maximum User Weight	100 kg	100 kg	125kg
Code	MFI-69397	MFI-70258	MFI-VNW01AYA

Sammons Prestons Wheelchair Foam Cushions

High-density Foam

Ideal for part- or full-time use

This cushion is perfect to fit into your wheelchair. It is lightweight and has ties on the back to allow for the cushion to be tied into the wheelchair for a secure fit and feel. Multiple sizes available to fit any wheelchair.

Key features

- High-density polyurethane foam and an extra strong poly cotton cover with ties
- Cover can be removed to launder
- Measures 16" and 18" deep
- Colour: Navy

Specifications

Pillow width	41 cm (16"), 46 cm (18")
Pillow depth	41 cm (16")
Pillow height	5 cm (2"), 7.6 cm (3")



	Pillow Size
PAT-081227461	41 x 41 x 5 cm (16 x 16 x 2")
PAT-081227479	41 x 41 x 7.6 cm (16 x 16 x 3")
PAT-081175082	46 x 41 x 5 cm (18 x 16 x 2")
PAT-081175090	46 x 41 x 7.6 cm (18 x 16 x 3")



Available colour






Rollators

Your guide to choosing the perfect rollator

Choosing a rollator with a frame at the correct height is essential to ensure wellbeing. Other factors to consider include where you will primarily use your rollator, the amount of support you will need, the number of wheels (three or four), portability and storage.

When taking measurements, consider the following:

- If the frame is too high, you may find it difficult to straighten your elbows and increase weight distribution through your arms.
- Your frame is too low if you are bent over with poor posture.
- Generally, the height of the handgrips should be at the level of the wrist bone when your elbows are slightly bent, at an angle of about **15° flexion**.
- When measuring the height of your walking frame, wear appropriate and supportive footwear.
- If your rollator includes a seat, consider the seat height. When seated comfortably, your feet should be flat on the ground.





Days 102 & 105 Rollator

Keeping you on the move

The 102 and 105 are our most popular rollators. With aluminium frames, combining lightweight design and robust engineering quality. Days 102 is compact with 6" wheels and Days 105 is standard with 8" wheels.

Key features

- Available in two sizes to accommodate a range of heights and builds
- Lockable loop brakes for security and convenience
- Adjustable brake tension to suit your needs
- Folds easily for storage and transport
- Ergonomic handgrips and adjustable height handles with five settings at 2.5-cm intervals
- ISO 11199-2:2005



	102	105
A	Handle height 77.5 - 87.5 cm (30½ - 34½")	Handle height 84 - 94 cm (33 - 37")
B	Seat height 50 cm (19¾")	Seat height 56.5 cm (22¼")
C	Backrest height 75 cm (29½")	Backrest height 80 cm (31½")
D	Wheel diameter 15 cm (6")	Wheel diameter 20 cm (8")
	Seat width 35 cm (13¾")	Seat width 32 cm (12¾")
	Product weight 6 kg (13.2 lbs)	Product weight 7.2 kg (15.8 lbs)



DAYS-105-BLUE
Available in 4 colours



DAYS-105-BLACK



DAYS-105-CHAMP




DAYS-102-RED
Available in Red

Colour	X-Small 102	Medium 105
Ruby Red	DAYS-102-RED	DAYS-105-RED
Blue		DAYS-105-BLUE
Black		DAYS-105-BLACK
Champagne		DAYS-105-CHAMP

<p>125kg 19½ st 275½ lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>2 YEAR Main frame warranty</p>
--	------------------------------------	---------------------------------------

Available colours



Rollators

Days Superlite

Quality and value


Deluxe super lightweight rollator with sleek and stylish Euro design. Large front wheels provides better stability and is ideal for outdoor walks. This lightweight, sturdy rollator collapses in one swift motion and stands for easy storage. Convenient push-button mechanism allows for quick height adjustment and hidden brake cable design provides added safety.

Key features

- Internal braking cable
- Lockable loop brakes for security and convenience
- Side folding for transport and storage
- Seat and storage bag
- Built in rest seat with canvas backrest.
- Walking stick or umbrella holder
- ISO 11199-2:2005



HCT-9137DR



A	Handle height 85 - 97 cm (33½ - 38¼")
B	Overall depth 70.5 cm (27¾")
C	Seat height 54 cm (21¼")
D	Backrest height 29.5 - 42 cm (11½ - 16½")
E	Front wheel diameter 25.4 cm (10") Rear wheel diameter 20.3 cm (8")
	Seat width 46 cm (18")
	Seat depth 25 cm (9¾")
	Product weight 6.91 kg (15¼ lbs)



Colour	
Ruby Red	HCT-9137DR
Blue	HCT-9137DB

<p>136 kg 21½ st 300 lbs Maximum user weight</p>	<p>1 YEAR Product warranty</p>	<p>2 YEAR Main frame warranty</p>
--	------------------------------------	---------------------------------------

Available colours



HCT-9137DB

Rollator in a Box

Standard and Heavy Duty

These standard and heavy duty rollators comes in a retail box with carry handle. Featuring 7½" wheels, it is one of the only rollators in the market with both height adjustable handles and seat height which makes it ideal for taller people. Easy to assemble with numbered snap in wheels and handbrakes.


Key features

- Easy to assemble - no tools required
- Retail carrying box with handle
- Height adjustable seat and handles
- Comes with bag
- Fold down backrest
- 7½" wheels
- ISO 11199-2:2005



KD1012

KD1014



	Standard KD1012 - 136kg	Heavy Duty KD1014 - 158 kg
A	Handle height 78.5 - 100 cm (30¾ - 39½")	Handle height 78.5 - 100 cm (30¾ - 39½")
B	Overall depth 68 - 80 cm (26¾ - 31½")	Overall depth 68 - 80 cm (26¾ - 31½")
C	Seat height 50 - 60 cm (19¾ - 23½")	Seat height 50 - 60 cm (19¾ - 23½")
D	Wheel diameter 19 cm (7½")	Wheel diameter 19 cm (7½")
	Seat width 35.6 cm (14")	Seat width 48.5 cm (19")
	Seat depth 30.5 cm (12")	Seat depth 30.5 cm (12")
	Overall width 62.2cm (24½")	Overall width 68.5cm (27")
	Product weight 7.4 kg (16½ lbs)	Product weight 7.7 kg (17 lbs)

Colour	Standard	Heavy Duty
Red	DAYS-KD1012-RED	DAYS-KD1014-RED
Blue	DAYS-KD1012-BLUE	DAYS-KD1014-BLUE

<p>Standard</p> <p>136 kg</p> <p>21½ st 299¾ lbs</p> <p>Maximum user weight</p>	<p>Heavy Duty</p> <p>158 kg</p> <p>24¾ st 348½ lbs</p> <p>Maximum user weight</p>	<p>1 YEAR</p> <p>Product warranty</p>
---	---	---------------------------------------

Available colours



KD1012

KD1014

Rollators

Days Seat Walker with Handbrakes Rollator

Ideal for everyday use

Height adjustable, lightweight, aluminium walker with lockable loop lever handbrakes and soft rubber hand grips. Features a curved padded backrest for comfort and a cushioned seat, which lifts to allow access to the handy shopping basket. Six inch wheels for a small turning circle allowing ease of use for both indoors and outdoors.

ARTG 199708

*Now ordered in Box Qty

Colour		Qty/Box
Red	MFI-MOBWAL70149-2	2
Blue	MFI-MOBWAL70493-2	2
Black	MFI-MOBWAL70168-2	2
Champagne	MFI-MOBWAL70425-2	2
Purple	MFI-MOBWAL70494-2	2



A	Handle height 83 - 95 cm (32¾ - 37½")
B	Overall depth 60 cm (26½")
C	Seat height 55 cm (21½")
D	Backrest height 82.5 cm (32½")
E	Wheel diameter 15 cm (6")
	Seat width 44 cm (17¼")
	Seat depth 30 cm (11½")
	Product weight 6.4 kg (14¼ lbs)

Available colour

Red

100 kg
 15¾ st 220½ lbs
 Maximum user weight

1 YEAR
 Product warranty

Days Seat Walker with Compression Brakes Rollator

Ideal for everyday use

A lightweight aluminium folding walker with compression (pushdown) brakes and soft rubber hand grips allowing easy cleaning to maintain hygiene, a feature not commonly found with budget walkers on the market. Colour red.

ARTG 199708

*Now ordered in Box Qty

Colour		Qty/Box
Red	MFI-MOBWAL70371-2	2
Rubber Tip	NOV-TA25-001	

Available colour


Red

100 kg
 15¾ st 220½ lbs
 Maximum user weight

1 YEAR
 Product warranty



A	Handle height 81 - 93 cm (31¾ - 36½")
B	Overall depth 58 cm (22¾")
C	Seat height 53 cm (20¾")
D	Backrest height 72.5 cm (28½")
E	Wheel diameter 15 cm (6")
	Seat width 30 cm (11½")
	Seat depth 30 cm (11½")
	Product weight 6.5 kg (14¼ lbs)



Days V4208 - Compact Seat Walker Rollator

Durable for everyday use

The Days V4208 Compact Seat Walker has a similar design to the Days V4206 Walkers, but is a more economical and compact model. This is a lightweight, stable walker that is easy to manoeuvre. It is easy to fold, making it easily transportable and features "Feather Touch" Handbrakes. It also has a large seat with extra padding, as well as a shopping compartment under the seat - making this a very versatile walker.

ISO 11199-2:2005

ARTG 199708

*Now ordered in Box Qty

Colour		Qty/Box
Red	MFI-V4208AUB-2	2
PVC Bag	MFI-SP-2000	

Available colour


Red

100 kg
 15¾ st 220½ lbs
 Maximum user weight

1 YEAR
 Product warranty



A	Handle height 73 - 84 cm (28¼ - 33")
B	Overall depth 60 cm (23½")
C	Seat height 48 cm (19")
D	Backrest height 75.5 cm (28½")
E	Wheel diameter 15 cm (6")
	Seat width 30 cm (11½")
	Seat depth 30 cm (11½")
	Product weight 7 kg (15½ lbs)



Days V4206 Walker Rollator

Durable for everyday use

The 4206 is our premium "A" frame seat walker. This walker features high quality 8" castors for superior outdoor use, while still being ideal for manoeuvring indoors. A large seat with extra padding and an extra large shopping bag under the seat make this a versatile walker. Available in a standard 550mm seat height or a low seat of 450mm.

ISO 11199-2:2005

ARTG 199708

*Now ordered in Box Qty

Colour	18" Low Seat	22" Standard Seat	Qty/Box
Red	MFI-V4206-18INCH-2	MFI-V4206-22INCH-2	2
PVC Bag	MFI-SP-2000	MFI-SP-2000	

Available colour

Red

100 kg
 15¾ st 220½ lbs
 Maximum user weight

1 YEAR
 Product warranty



A	Handle height V4606: 70 - 79 cm (27½ - 31") V4206: 80 - 91 cm (31½ - 36")
B	Overall depth V4606: 64 cm (25¼") V4206: 68 cm (26¾")
C	Seat height V4606: 46 cm (18") V4206: 56 cm (22")
D	Backrest height V4606: 72.5 cm (28½") V4206: 82.8 cm (32½")
E	Wheel diameter 20 cm (8")
	Seat width 34.5cm (13½")
	Padded seat 35 mm x 33 mm x 38 mm
	Product weight 7 kg (15½ lbs)

Rollators

Bariatric Rollators

Built for strength

The 111 Bariatric Rollator, constructed of high-quality steel, was designed to help users maintain balance and prevent falls while navigating indoor and outdoor surfaces. The rollator is height-adjustable with a padded backrest to prevent the user from tilting backward when seated. Complete with a storage basket for personal belongings, this rollator folds for easy storage.

Key features

- High-quality steel frame with height-adjustable handgrips
- Lockable loop brakes with adjustable tension
- Extra-wide padded seat and backrest
- ISO 11199-2:2005
- ARTG 174476



*Now ordered in Box Qty

Specifications	Low Mack	Supa Mack	High Mack*
Weight of Walker	9.4kg	9.4kg	9.4kg
Floor to Seat Height	490mm	550mm	580mm
Seat Wide	500mm	610mm	500mm
Seat Depth	297mm	350mm	297mm
Handle Height Range**	750mm - 860mm	750 - 900mm	870 - 980mm
Overall Width (footprint)	700mm	710mm	700mm
Overall Depth (footprint)	640mm	650mm	640mm
Wheels - Front	200mm (8")	200mm (8")	200mm (8")
Wheels - Rear	200mm (8")	200mm (8")	200mm (8")
Maximum User Weight	150 kg	225 kg	150 kg
Code	NOV-AC34L	NOV-MOBWAL70116	NOV-AC34H-2

Available colour




Heavy-duty Rollator 220kg

Built for strength

Heavy-duty Rollators are manufactured with gauge steel tubing and offer wider frames and rest seats for users who need more space.

Key features

- Lockable loop brakes for security and convenience
- Under-seat basket for functional storage space
- High user weights
- Folds easily for transport and storage
- ISO 11199-2:2005



A	Handle height 83 - 94 cm (32¾ - 37")
B	Overall depth 66 cm (26")
C	Overall width 67 cm (26½")
D	Seat height 55 cm (21½")
E	Wheel diameter 20 cm (8")
	Seat width 49 cm (19¼")
	Seat depth 35 cm (13¾")
	Product weight 13 kg (28.7 lbs)

Colour	
Red	DAYS-HD-RED
Blue	DAYS-HD-BLUE



Available colours



Rollators


Breeze Indoor Rollator

A smooth journey from room to room

Breeze Indoor Rollators offer a stylish alternative to a household trolley. Manufactured from painted aluminium tubing with four solid tyre wheels, each rollator includes a storage bag, locking brakes and height-adjustable handles.

Key features

- One-handed braking with brake lock
- Easily adjustable handle height
- Durable plastic tray and carry bag with zip pocket
- Fold button for easy storage



A	Handle height 83.5 - 94.5 cm (32¾ - 37¼")
B	Overall depth 66 cm (30")
C	Overall width 54 cm (21¼")
D	Wheel diameter 17.75 cm (7")
	Tray depth 29 cm (11½")
	Tray width 39 cm (15½")
	Tray height 68 cm (26¾")
	Tray max load 5 kg (11 lbs)
	Bag max load 5 kg (11 lbs)
	Product weight 5.5 kg (12 lbs)



Brake lock



Easily adjust handle height

Colour	
White	PAT-091558360

Available colour



Adjustable Height Trolley Walker

Comfort when you need it


The adjustable trolley walker is supplied with a cream painted steel frame. Clip on trays can easily be removed for cleaning. Adjustable height easily accessible loop lockable brakes. As the trolley has brakes, this Trolley can also be a walking aid.

Key features

- Cream painted steel frame
- Handles are adjustable in height
- Plastic clip on trays, easily removable for cleaning
- Easily accessible loop lockable brakes
- Can be used with or without trays or can be used as a table
- ISO 11199-2:2005
- ARTG 199708

Specifications

Overall width	58.5 cm (23")
Overall depth	56 cm (22")
Handle height range	75.5 - 90.5 cm (29¾ - 35½")



A	Handle height 75.5 - 90.5 cm (29¾ - 35½")
B	Overall depth 56 cm (22")
C	Overall width 58.5 cm (23")
D	Wheel diameter 15 cm (6")
	Product weight 8.25 kg (18.2 lbs)



*Now ordered in Box Qty

Colour		Qty/Box
Beige	HLC-109-2	2
White	HLC-109NS-2	2
Tray	HLC-PL244	

Available colours



Tri-wheel Walkers


Aluminium Lightweight Tri-wheel Walker

Light and easy to manoeuvre

Aluminium Lightweight Tri-wheel Walkers feature a lightweight aluminium frame that is easy to manoeuvre, making them ideal for use indoors or outdoors.

Key features

- Lockable loop brakes for security and convenience
- Available with or without a basket
- Includes durable bag to transport and conceal your belongings
- ISO 11199-2:2005



A	Handle height 83 - 93 cm (32¾ - 36½")
B	Overall depth 61 cm (24")
C	Overall width 68.5 cm (27")
D	Wheel diameter 19 cm (7½")
	Product weight 7.7 kg (17 lbs)




Optional basket and tray



*Now ordered in Box Qty

Colour	With bag only	Qty/Box
Ruby Red	DAYS-TRI-RED-2	2
Blue	DAYS-TRI-BLUE-2	2
Accessories		
Basket and Tray	DAYS-091555242	

Available colours




2-in-1 Rollator and Transit Chair

Comfort when you need it

The 2-in-1 Rollator and Transit Chair offers the comfort you need for various uses. Each rollator includes a padded seat and backrest, an under-seat bag for storage and flip-up footrests to provide extra support.

Key features

- Suitable for use as a rollator or transit chair
- Lockable loop brakes for security and convenience
- Folds easily for transport and storage
- ISO 11199-2:2005



A	Handle height 80 - 93 cm (31¼ - 36½")
B	Pushing handle height 90.5 cm (35½")
C	Overall depth 74 cm (29¼")
D	Overall width 58.5 cm (23")
E	Seat height 56 cm (22")
F	Wheel diameter 10 cm (8")
	Seat width 46 cm (18¼")
	Seat depth 33 cm (13")
	Product weight 8.8 kg (19.4 lbs)



Colour	
Red	DAYS-2IN1-RED
Silver	DAYS-2IN1-SIL

Available colours



Standing Walkers


Gutter Walker

Sturdy, robust and economical

Designed to keep the user in an upright position for gait training. Folds for storage. Forearm support and hand are fully adjustable.

Specifications

Overall height	110 - 134 cm (43¼ - 52¾")
Overall forearm height	99 - 122 cm (39 - 48")
Width between wheels	58 cm (22¾")
Overall width	63 cm (24¾")
Overall depth	90 cm (35½")
Wheels	19 cm (7½")
Weight capacity	150 kg (330½ lbs)



A	Forearm height 98 - 122 cm (38½ - 48")
B	Overall height 110 - 134 cm (43¼ - 52¾")
C	Overall depth 90 cm (35½")
D	Overall width 63 cm (24¾")
E	Width between wheels 58 cm (22¾")
F	Wheel diameter 19 cm (7½")
	Product weight 11.5 kg (25.4 lbs)

Colour	
Grey	DAYS-203B

Available colour

Grey	150 kg 23½ st 330½ lbs Maximum user weight	1 YEAR Product warranty
------	--	----------------------------



Walking Frames

Go further with Days


All Days Walking Frames are manufactured from lightweight anodised aluminium that is durable and easy to wipe clean. The plastic handgrips have an ergonomic indented shape for greater comfort. Non-slip rubber ferrules are non-marking and provide excellent grip and stability.



PAT-AA8389	Walking Frame Caddy
------------	---------------------

Walking Frames

Ball Walker Walking Frame



125 kg
19 3/4 st 275 1/2 lbs
Maximum user weight

1 YEAR
Product warranty

Additional Options

JAN-DAYBALLWALKGL - Pair of Glides
JAN-DAYBALLWALKWH - Pair of Wheels



	A - Height	B - Overall width	C - Depth	Product weight
JAN-DAY BALL WALKER	87 - 92 cm (34 1/4 - 36 1/4")	44 - 64 cm (17 1/4 - 25")	35 - 44 cm (13 3/4 - 17 1/4")	2.0 kg (4 1/4 lbs)

ARTG 174419

Side Fold Walking Frame



125 kg
19 3/4 st 275 1/2 lbs
Maximum user weight

1 YEAR
Product warranty

Additional Options

JAN-JMC-3725 - Pair of Glides
JAN-JMC-3700 - Pair of Wheels



Without wheels	Qty/Box	A - Height	B - Overall width	C - Depth	D - Handle to handle	Product weight
JAN-JMC-233EL-2	2	81 - 91 cm (31 3/4 - 35 3/4")	55 cm (21 1/2")	45 cm (17 3/4")	43 cm (17")	2.7 kg (5 3/4 lbs)

ARTG 174419



Crutches

All Days Crutches are made from lightweight aluminium with a tough anodised finish, making them incredibly durable and easy to use. They are designed for minimum maintenance and only require an occasional wipe-down with a soft cloth. Disinfectant solution can be used if needed.



Crutches

How to select crutches

Factors include handle options, adjustment levels and standard features

Handle options



Standard

Standard plastic handgrips are easily wiped clean to reduce the risk of cross-contamination.



Ergonomic

Ergonomic handles were designed and developed with physiotherapists and healthcare professionals. The unique handle distributes pressure over a greater area of the user's hand, improving comfort and safety.



Comfy grip

Comfy grip handles are made from soft-touch plastic to provide a more comfortable and cushioned grip — ideal for long-term use.

Adjustment levels



Single adjustable crutches

Single adjustable crutches enable you to adjust the handle height, but the distance between the handle and arm cuff is fixed.



Double adjustable crutches

Double adjustable crutches enable you to adjust the handle height and the distance between the handle and arm cuff.



Permanent crutches

Permanent crutches are not adjustable and must be cut to size to suit each user. They are ideal for long-term use.

Standard features



Closed fixed armband

The closed fixed armband enables you to retain the crutch on your arm even if you release your grip on the handle. This prevents the crutch from dropping to the floor accidentally.



Non-marking ferrules


Non-marking ferrules are slip-resistant for added safety. The 22-mm ferrules are fitted for adult crutches, while 16-mm ferrules are fitted for children's crutches.



Push-pin adjustable

Adjustment is easy using the simple push-pin operation.

Forearm Crutches*



		A - Handle to armband	B - Handle to ground
Single adjustable			
JAN-125A	Medium	23 cm (9")	89 - 120 cm (35 - 47 1/4")
JAN-125AT	Large	23 cm (9")	89 - 137 cm (35 - 54")

ISO 11334.1.2014
ARTG 174418


Ergonomic Grip Forearm Crutches*



		A - Handle to armband	B - Handle to ground
Double adjustable			
JAN-124A	Medium	20-27.5 cm (8-10 1/2")	89 - 120cm (35 - 47 1/4")
JAN-124AT	Large	20-27.5 cm (8-10 1/2")	89 - 137 cm (35 - 54")

ISO 11334.1.2014
ARTG 174418

Aluminium Underarm Crutches*



		A - Armband to ground
Double adjustable		
DAYS-80536	Small	94 - 114 cm (37 - 44 3/4")
DAYS-80535	Medium	114 - 135 cm (44 3/4 - 53")
DAYS-80534	Large	135 - 155 cm (53 - 61")

ISO 11334.1.2007
ARTG 277315

Additional Options



Crutch Tips

Spare rubber tips to suit all aluminium underarm crutches

	Colour	Qty/Size
JAN-PL155	Grey	Pair / 22cm



Crutch Pads

Spare rubber pads to suit all aluminium underarm crutches

	Colour	Qty
JAN-JMC2014	Grey	Pair



Crutch Grips

Spare rubber grips to suit all aluminium underarm crutches

	Colour	Qty
JAN-JMC2015	Grey	Pair

*Crutches are sold as pairs.



Walking Sticks

Select the walking stick that is right for you


Days Walking Sticks are available in various styles. Use the guide on the following page to help you select the right walking stick for you.



How to select walking sticks

Factors include stick height, left- or right-handed and type

Walking stick height




Using a walking stick at the correct height is essential. When measuring, the handle should rest on the floor while your arm is relaxed. Then, adjust or trim the stick to the **crease of the wrist**.

When using your walking stick, you should maintain an upright posture with the elbow slightly bent. Your body weight will be supported when you push down on the stick while walking.


If you need two walking sticks, the length should be slightly longer as they will be used ahead of the body when you are walking.

Left or right hand



You should hold the stick in hand on the **opposite side to the one that needs support** if you are using a walking stick to compensate for an injury to your leg or joint on one side. If you use the stick for balance, you can use it with either hand.

Types of walking sticks



Adjustable

Adjustable walking sticks can be modified to different heights using simple spring-loaded catches and telescopic shafts.



Tripod and tetrapod

Tripod and tetrapod walking sticks have a broad footprint to provide additional stability. These are best suited for heavier users who require a more stable base.




Folding walking sticks

Folding walking sticks are lightweight metal sticks with sectioned shafts that enable you to fold them for storage. For example, you may transport them within a handbag. Strong elastic runs inside the shaft to ensure that the stick remains stable in its open position.

Walking Sticks

T-Handle Walking Stick



113 kg
17³/₄ st 249 lbs
Maximum user weight

1 YEAR
Product warranty

		A-Handle to ground
DAYS-20052BR	Bronze	78 - 101 cm (30 ³ / ₄ - 39 ³ / ₄ ")

ARTG 277315

Swan Handle Walking Stick



113 kg
17³/₄ st 249 lbs
Maximum user weight

1 YEAR
Product warranty

		A-Handle to ground
DAYS-20042SI	Silver	81 - 103.5 cm (31 ³ / ₄ - 40 ³ / ₄ ")

ARTG 277315

Anatomical Walking Stick



113 kg
17³/₄ st 249 lbs
Maximum user weight

1 YEAR
Product warranty

	Hand	A-Handle to ground
DAYS-20051L	Left	79 - 101 cm (31 - 39 ³ / ₄ ")
DAYS-20051R	Right	79 - 101 cm (31 - 39 ³ / ₄ ")

ARTG 277315

Folding Walking Stick



113 kg
17³/₄ st 249 lbs
Maximum user weight

1 YEAR
Product warranty

		A-Handle to ground
DAYS-20061LBK	Black	79 - 101 cm (31 - 39 ³ / ₄ ")

ARTG 277315

Walking Stick Wrist Straps



NOV-WRISTSTRAPS	10 Pack	21.5cm
-----------------	---------	--------

Black Rubber Ferrules



PAT-091168178	Each	13 mm
PAT-AA8200CW	Pair	19mm

Walking Sticks

Crook Handle Walking Stick



	A-Handle to ground
DAYS-20720BR	73 - 96 cm (28 ³ / ₄ - 37 ³ / ₄ "

ARTG 277315


Quadruped Walking Stick



	Base size	A-Handle to ground
JAN-JMC-2660	Small	69 - 94 cm (27 - 37")
JAN-JMC-2640	Wide	74 - 100 cm (29 - 39")

ARTG 174419 ISO 1134-4:1999

Folding Patterned Walking Stick



	Pattern	A-Handle to ground
DAYS-P193	Floral	84 - 94 cm (33 - 37")
DAYS-P081	Blossom	84 - 94 cm (33 - 37")
DAYS-P678	Paisley	84 - 94 cm (33 - 37")
DAYS-P554	Leopard	84 - 94 cm (33 - 37")
DAYS-P198	Pastel	84 - 94 cm (33 - 37")
DAYS-PMIX	Mixed x25 pk	84 - 94 cm (33 - 37")

ARTG 277315

Tripod Seat Stick



	A-Handle to ground	B-Seat to ground
NOV-MOBWAL68073	71.5 cm (28")	54 cm (21 ¹ / ₄ "



How to order

AUSTRALIA

Performance Health ANZ
P: 1300 473 422
E: sales@performhealth.com.au
W: performhealth.com.au
Unit 3, 3 Basalt Road
Pemulwuy, NSW 2145

UK CUSTOMER SERVICE

P: 03448 730 035
F: 03448 730 100
E: UKCustomerSupport@performancehealth.com
08:30-16:00 (CET)

UK OFFICE

Performance Health International Ltd.
Nunn Brook Road,
Huthwaite,
Sutton-in-Ashfield,
Nottinghamshire,
NG17 2HU

INTERNATIONAL CUSTOMER SERVICE

P: +44 1623 448 706
F: +44 1623 448 784
E: international@performancehealth.com

Worldwide locations

US CORPORATE HEAD OFFICE

28100 Torch Parkway, Suite 700
Warrenville, IL 60555
P: 630-393-6000
F: 630-393-7600

FRANCE

13 Rue André Pingat
CS 10045
51724 Reims Cedex
France
P: +33 3 10 00 79 30
E: fr-contact@performancehealth.com

GREECE

208 Kifissias Avenue
15231 Halandri
Athens, Greece
P: +30 210 6712 266

Days® is a registered trademark of Performance Health and/or its subsidiaries and may be registered in the United States and other countries.
©2022 All rights reserved. Unauthorized use is strictly prohibited.

DAYS22CAT

